

AL-LPO 1

Bottle Strain Viewer-Light Polarizer

Inspects and catches defects of Bottle Preforms before allowing abundance and wasteful units entering into Blow Molding process.
Irregularities and Residual Stress/Strains are shown clearly



A standard model holds 10 piece of Preforms for 16 ounce or 600 ml Bottles. Other configurations and fixture sizes are specified at time of ordering
More or less fixtures in a single inspection pass depending on size of Preforms

at **AdvanceLab**, we pride ourself with producing instruments with **Highest Quality Standards**. Fully featured designs such as :

Continuous Motorized Automatic run & **Jog** to specific position

Manual Hand turning, with 90° indexing stops

* (other sizes and multiple level larger systems and custom engineered versions are available, please ask our sales for details) *

Technical Details

| Motorized Automatic | Accommodation | Manual Hand Turning Indexing |
|---|---|---------------------------------|
| continous motor rotation & Jog Indexing | standard model holds 10 preforms in single pass | full 360° or 90° indexing stops |
| 110/ 220 VAC ; 5 amps ; 50/60 Hz | | |

Shipping Weight : 44 lbs (20 kgs) • Packing : 19" x 24" x 15" (48 x 61 x 38cm) • Net Weight : 40 lbs (18 kgs)